

ETHANOL

Enzymatic UV


Summary:

The determination of ethanol belongs to the most frequent analyses in the forensic and toxicological laboratory. It serves for the diagnosis of intoxications and poisonings particularly for emergency room patients.

Method:

Enzymatic UV test with alcohol dehydrogenase (ADH)

Principle:



In the presence of NAD Ethanol is converted by the alcohol dehydrogenase. The measured absorbance of the produced NADH is proportional to the ethanol concentration in the sample

Reagents:

Components and Concentrations

N.B. Concentrations are those in the final test mixture

R1: Buffer pH 9.0 240 mmol/l

Stabilizers and preservatives

R2: Buffer Ph 6.6 8 mmol/l

NAD ≥2 mmol/l

Alcohol dehydrogenase ≥40 kU/l

Storage Instructions And Reagent Stability

The reagents are stable up to the end of the indicated month of expiry, if stored at 2-8 °C and Contaminations is avoided. Do not freeze the reagents! The standards are stable up to the end of the indicated month of expiry, if stored at 15 – 25 °C.

Warnings And Precautions

1. The reagents contain sodium azide (0.95 g/l) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
2. Take the necessary precautions for the use of laboratory reagents

Waste Management

Please refer to local legal requirements.

Reagent Preparation

The reagents and the standards are ready to use.

Materials Required But Not Provided

NaCl solution (0,9 %).

General laboratory equipment

Assay Procedure

Application sheets for automated systems are available on request.

Wavelength 376 nm(360-380)

Optical path 1 cm

Temperature 37 °C

Measurement Against reagent blank

The observance of exact measuring times and absolute equal treatment of all samples, standards and controls must be respected

Reagent blank Sample /Standart

Sample / Standart – 10 µl

Dist. Water 10 µl –

Reagent I 1000 µl 1000 µl

Mix and incubate 5 min. at 37 °C Read absorbance A1 then add

Reagent II 1000 µl 1000 µl

Mix and incubate 5 min. at 37 °C Read absorbance A2 immediately.

A= (A2 – A1) Sample /Standart

Specimen

Serum and plasma (heparin and EDTA)

The stability in serum and plasma is 2 weeks at 20-25 °C, 6 months at 4-8 °C. Samples must be stored tightly closed. Don't use alcohol or volatile disinfectants during ethanol measurement. Discard contaminated specimens!

Calculations

One point or multi-point calibration.

One point calibration:

With standart 1.0 g/l ;

A Sample

x C Standard (g/l)

A Std

Multi-Point Calibration

The ethanol concentration of unknown samples is derived from a calibration curve using a linear Algorithm. The calibration curve is obtained with four standards at different levels (see order information) and NaCl solution (0,9 %) for determination of the zero value.

Conversion Factor

Ethanol (g/l) *21.7 = Ethanol (mmol/l)

Ethanol (g/l)* 0.8 = Ethanol %

Performance Characteristics

| Intra assay n=20 | Mean g/l | SD g/l | CV % |
|---------------------|-------------|-----------|---------|
| Sample 1 | 0.51 | 0.01 | 1.67 |
| Sample 2 | 0.98 | 0.02 | 1.95 |
| Sample 3 | 1.99 | 0.01 | 0.66 |
| Inter assay n=20 | Mean g/l | SD g/l | CV % |
| Sample 1 | 0.51 | 0.02 | 3.36 |
| Sample 2 | 1.01 | 0.02 | 2.05 |
| Sample 3 | 1.99 | 0.03 | 1.66 |

Measuring Range

The test has been developed to determine ethanol concentrations up to 3.5 g/l. When values exceed this range samples should be diluted 1+1 with NaCl solution (0,9 %) and the result multiplied by

Specificity interferences

No interference was observed by ascorbic acid up to 30 mg/dl, bilirubin up to 60 mg/dl, lipemia up to 2000 mg/dl triglycerides, hemoglobin up to 1000 mg/dl, creatinine up to 250 mg/dl glucose up to, 2000 mg/dl, urea up to 2000 mg/dl and LDH up to 2000 U/l.

Sensitivity / Limit of Detection

The lower limit of detection is 0.1 g/l.

Precision (at 37 °C)

Method Comparison

A comparison between BIOANALYTIC Ethanol (y) and a commercially available assay (x) using 30 samples gave following results : y= 1.0 * + 0.1 g/l ; r = 0.999.

Reference Range

Result < 10 mg/dL Negative result
 10 – 30 mg/dL 2.17 – 6.5 mmol/L Without obvious signs
 30 – 120 mg/dL 6.5 – 26.0 mmol/L Slowed reflexes, diminution of attention, judgment and control
 120 – 250 mg/dL 26.0 – 54.3 mmol/L Reduced visual acuity and increased reaction time
 250 – 350 mg/dL 54.3 – 76.0 mmol/L Muscular incoordination, decreased response to stimuli
 > 3.5 g/L > 76.0 mmol/L Impairment of circulation and respiration, possible death



ETHANOL

Enzymatic UV

Quality Control:

Control Serum: BIOANALYTIC ETH CON.

SET 5 x 1 ml #B10827

The control intervals and limits must be adapted to the individual laboratory and country-specific requirements. Values obtained should fall within established limits. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

Calibration:

Standardization: This Ethanol method was calibrated against an international standard defined for Ethanol.

S1: BIOANALYTIC ETH CAL. SET 5 x 1 ml #B11914

Calibration frequency:

Calibration is recommended after lot change as required following quality control procedures

Bibliography

1. Borque L, Rus A, Dubois H. Automated determination of streptolysin O antibodies by turbidimetric latex immunoassay method. J Clin Immunoassay 1992; 15: 182-6.
2. Thomas L. Clinical Laboratory Diagnostics I ed Frankfurt: TH-books Verlagsgesellschaft, 1998. p.1168 – 1170 .
3. Immunology and Serology in Laboratory Medicine, 2nd edition. Turgeon ML. Mosby, 1996.








Order information (Cat No.) :

| | | | | | |
|--------|--------|--------|--------|--------|--------|
| CR535 | B24146 | B28145 | B31146 | B34145 | B36146 |
| B21145 | B25145 | B28146 | B32145 | B34146 | B37145 |
| B21146 | B25146 | B30145 | B32146 | B35145 | B37146 |
| B22145 | B27145 | B30146 | B33145 | B35146 | B80145 |
| B24145 | B27146 | B31145 | B33146 | B36145 | |

Manufacturer

Diaclinica Diagnostik Kimya.San.Tic.Ltd.Şti
Adress : İkitelli O.S.B Mutsan San.Sit. M4 Blok No:19Başakşehir/İSTANBUL
Tel:+90(212) 549 33 88- Fax:+90 (212) 549 55 50
Web :www.diaclinica.com

SYMBOLS

| | |
|---|----------------------------------|
|  | for in vitro diagnostic use only |
|  | lot of manufacturing |
|  | code number |
|  | storage at temperature interval |
|  | expiration date (year/month) |
|  | warning, read enclosed documents |
|  | Read the directions |



ISO 9001:2015
ISO 13485:2016

